

# 1000 Shades of White

... all in 1 Syringe!

Paste

New:  
Flowable



Tokuyama Dental

omniCHROMA *Flow*

The Future of Composites and Flowables:  
Colour through Light

From the Technology Pioneer.



Tokuyama Dental  
Dental High Tech from Japan

# Unique worldwide: Only a single composite/flowable for the entire range of tooth shades

After OMNICHROMA for composites, TOKUYAMA DENTAL has once again made the decisive innovative leap: this time with OMNICHROMA FLOW in the field of flowables. For the second time, the long sought-after chameleon effect has been achieved successfully culminating in natural perfection. The reason is as simple as it is spectacular: OMNICHROMA and OMNICHROMA FLOW as well as the chameleon are coloured through light.

Tokuyama Dental

**omniCHROMA**

Tokuyama Dental

**omniCHROMA Flow**



In conventional composites, a limited number of tooth shades, for example, from A1 to D4 according to the VITA shade system, are reproduced using colour pigments. Due to the “*Smart Chromatic Technology*”, OMNICHROMA and OMNICHROMA FLOW completely dispense with colour pigments and instead use the natural principle of structural colour – colour that becomes visible when light strikes special structures, the same as with chameleons.

- In 2015, scientists at the University of Geneva discovered that the chameleon has an omentum of nanocrystals in its skin that selectively reflects certain wavelengths of light.

- OMNICHROMA and OMNICHROMA FLOW consist of a homogeneous “pearl structure”, which makes the reflection of a precisely defined light wave range possible.
- The targeted refraction of light creates structural colour in the yellow-red range and also reflects the surrounding real tooth colour.
- This was achieved by further developing TOKUYAMA DENTAL’s patented “*Sub-Micro-Pearl-Technology*” with spherical fillers obtained according to the “cultured pearl principle”.

## Spectrum of indications

- Direct anterior and posterior restorations
- Directly bonded composite veneers
- Diastema closure
- Repair of ceramic/composite



### Numerous advantages: OMNICHROMA and OMNICHROMA FLOW


- Eliminates the need for shade determination
- Simplifies stocking
- Eliminates the need for special colours
- Reduction in expiring material
- Permanent availability of the right shade

# Uniquely aesthetic: Outstanding colour adaptation

Smart Chromatic Technology: precisely 0.00026 millimeter sized beads and the translucency of the cured composite ensure that the right tooth shade is always created through light.


Conventional composites

**THE COLOUR PIGMENTS  
OF THE COMPOSITE  
(FOR EXAMPLE, A3)  
ARE REFLECTED.**



OMNICHROMA

**THE REAL TOOTH SHADE  
OF THE CAVITY WALLS  
OR CAVITY BASE IS  
REFLECTED AS IS  
THE RED-YELLOW  
STRUCTURAL COLOUR  
OF THE HOMOGENEOUS  
PEARL STRUCTURE.**



**The innovation  
breakthrough from  
TOKUYAMA DENTAL's  
research**

The phenomenon of colour adaptation in chameleons is the same as with OMNICHROMA: the light strikes a very specific structure in the sub-micro or nanometer range.

The decisive question was: which structure, in other words, which size and shape do the fillers need to have for the light to reflect real tooth shades optimally?

In 2018, TOKUYAMA DENTAL's research found the answer: the "Smart Chromatic Technology". When light strikes small spherical fillers at exactly  $0.26 \mu\text{m}$  the refraction and diffraction of the light generates the ideal red-yellow colour effect that is necessary to faithfully imitate the genuine tooth shade.

In contrast to conventional composites, OMNICHROMA is intended to achieve an extreme colour change after curing and is also very practical when modelling.



Typical for OMNICHROMA:  
Opaque beforehand ...  
beautifully semi-translucent  
afterwards

SM

CHRO  
TECHN

As many  
white a  
are p



Johanna's white



Leo's white



Betty's white



# ART MATIC OLOGY

shades of  
s there  
eople.

The key technology for OMNICHROMA comes from the TOKUYAMA DENTAL Research Center in Japan.



The TOKUYAMA DENTAL Corporation, based in Japan, has been developing

innovative dental solutions for over 40 years and is one of the leading manufacturers of products for conservative and prosthetic dentistry. In 2020, this innovative strength was recognized internationally. OMNICHROMA and TOKUYAMA DENTAL received awards as the most innovative material and most innovative company from the renowned Dental Advisor.



Sophie's white




Anna's white



# The secret is our way of manufacturing


## Conventional process

### PRODUCTION OF FILLERS BY A GRINDING PROCESS




## Sub-Micro-Pearl-Technology

### “GROWING” FILLERS BY THE SOL-GEL METHOD




Glass materials are ground until the individual filler particles are within a desired size range. However, an exact control is not possible in this way, so that the filler particles differ significantly in shape and size. This becomes particularly clear as soon as the fillers in question are observed under the scanning electron microscope.


TOKUYAMA DENTAL produces OMNICHROMA fillers based on its own patented “*Sub-Micro-Pearl-Technology*”. In this process, the Sol-Gel method is used to progressively coat spherical fillers in an organic solution. After several weeks, the fillers have “grown” evenly in a spherical shape and are exactly 0.26 µm in size. In this optimal size, the desired colour adaptation effect is achieved precisely in combination with other outstanding physical properties.



Venus One, Kulzer  
(1 µm; 20,000x magnification)



Tetric EvoCeram, Ivoclar Vivadent  
(1 µm; 20,000x magnification)




OMNICHROMA, TOKUYAMA DENTAL  
(1µm; 20,000x magnification)



Filtek Supreme XTE Flow, 3M Espe  
(1 µm; 20,000x magnification)



G-aenial Universal Flow, GC  
(1 µm; 20,000x magnification)



OMNICHROMA FLOW, TOKUYAMA DENTAL  
(1µm; 5,000x magnification)





## Clinical studies confirm the excellent colour adaptation properties of OMNICHROMA (CAP-V)

The University of Texas study analysed the Visual Colour Adjustment Potential (CAP-V) of various composite materials through visual evaluation.

Of the five composites tested, OMNICHROMA from TOKUYAMA DENTAL demonstrated the best shade adaptation effect.

The evaluation of shade differences in Class I restorations compared to the surrounding artificial tooth substance was lowest for OMNICHROMA, which represented the best correlation between OMNICHROMA and the A1-D4 prosthetic teeth.



Center: Denture tooth with Class I cavity


Circle: A1-D4 denture teeth restored with OMNICHROMA (dual specimen)

Evaluation of Color Adjustment Potential of Resin Composites Pereira-Sanchez N., Paravina R.D., et al. (University of Texas)

## Clinical studies confirm the excellent colour adaptation properties of OMNICHROMA (CAP-I)

Another study by the University of Texas also analysed the Instrumental Colour Adjustment Potential (CAP-I) of various composite materials by instrumental evaluation.

OMNICHROMA again confirmed the results of the previous study (CAP-V) in the instrumental determination of the colour adaptation effect. Thus OMNICHROMA offers both visually and instrumentally a very broad chameleon effect across the entire VITA shade range.



# OMNICHROMA:

## International study results

### Effect of filler size on the colour adaptation effect

To demonstrate that the size of the fillers significantly relates to the structural colour produced, Tokyo Medical and Dental University examined three composites with fillers of different sizes for their respective shade adaptation effects. Cavities on eleven different human teeth (shades A2, A3, A4, B2, B3, B4, C2, C3, C4, D2, and D4) were filled with OMNICHROMA (260 nm), an experimental nano-composite (100 nm) and a conventional composite with ground fillers. The respective colour adaptation was then measured using a spectrophotometer and compared with the measured colour before filling (Delta E<sub>00</sub>).


	ECM	R1	R2
.....			
FROM INCISAL	1.40±0.37 <sup>A</sup>	2.09±73 <sup>B</sup>	2.02±0.68 <sup>B</sup>
.....			
FROM CERVICAL	1.45±0.30 <sup>a</sup>	2.07±0.23 <sup>b</sup>	2.05±0.21 <sup>b</sup>

The samples demonstrated significant differences with respect to Delta E<sub>00</sub> ( $p < 0.05$ ), depending on the material differences. OMNICHROMA exhibited significantly lower Delta E<sub>00</sub>-values than R1 (ground fillers) and R2 (nano fillers) from both viewed from incisal and cervical, suggesting that OMNICHROMA has a superior ability to mimic the original shade of the human tooth. This also indicates that the particle size of the filler has a significant influence on OMNICHROMA's colour adaptation capability.

### Not only superior in colour adaptation, but also superior under load

The Oregon Health & Science University investigated numerous material properties of OMNICHROMA in comparison to conventional composite materials. Here, too, OMNICHROMA's superior class was demonstrated, as the special filler structure and filler composition result in an outstanding load-bearing capacity of the material. Among other things, this is reflected by excellent compressive strength.

OMNICHROMA FLOW also proves to be extremely resilient in tests. If paste-like composites are generally considered to be more resilient due to their higher filler content, the flowable OMNICHROMA FLOW with almost 400MPa compressive strength does not have to shy away from comparison with the supposedly stronger composites. OMNICHROMA FLOW thus opens up the entire range of indications to the user without any compromises.





## **Not only superior in colour adaptation, but also superior colour retention**

At Okayama University, a series of tests investigated the colour stability of OMNICHROMA. Four artificial teeth in the shades A1, A2, A3 and A3.5 were filled accordingly and stored in a water bath. The colour adaptation of OMNICHROMA was measured immediately after one day, after one week and after three months. The result was as follows: "The study demonstrated that OMNICHROMA can cover a wide range of cavity colours with just a single shade. This study also demonstrated that the composite material based on the structural shade system has the potential to enable aesthetic restorations without shade selection or layering."

A further in-house study (TOKUYAMA DENTAL R&D) also showed that OMNICHROMA still adapts to the changed real tooth colour even after bleaching.



# The surface is decisive


Which is easier to polish?  
Which reflects light more uniformly?  
Which shines more beautifully?

# Brilliant results


With its surfaces smooth like mirrors, TOKUYAMA DENTAL's Sub-Micro-Pearl-Technology provides a fast and long-lasting gloss.

In the case of typical ground fillers, the light strikes extremely irregular surfaces that scatter the light diffusely and appear correspondingly matt or require a long and complex polish.

The law of reflection "angle of incidence equals angle of reflection" is the basic prerequisite for the gloss effect and only works with very smooth surfaces: with mirrors, natural teeth and OMNICHROMA.




CONVENTIONAL NANO-HYBRID FILLERS  
ROUGH, IRREGULAR SURFACES REFLECT LIGHT MORE DIFFUSELY: LESS GLOSS





OMNICHROMA WITH SUB-MICRO-PEARL-TECHNOLOGY  
UNIFORM REFLECTION OF LIGHT AS WITH A MIRROR: NATURAL HIGH GLOSS

**Superior polishability:**  
**89 % gloss already after 30 sec.**



**Excellent abrasion properties**

OMNICHROMA is a particularly abrasion-resistant composite, yet at the same time gentle on the antagonist due to its unique filler structure, which only offers a small surface area for abrasion.



Filtek Supreme XTE,  
3M Espe

before and after the  
abrasion test  
(50.000 cycles)



OMNICHROMA,  
TOKUYAMA DENTAL

before and after the  
abrasion test  
(50.000 cycles)

TOKUYAMA DENTAL R&D

# OMNICHROMA: A composite with unique properties



- Due to the perfect chameleon effect, time-consuming shade determinations are a thing of the past
- Always reorder, store and document only 1 product at a time
- Highest level of polishability: glossy like a mirror in a flash due to spherical fillers

- Creation of a structural colour through precise light refraction of the homogeneous pearl structure as well as reflection of the surrounding tooth shade: for a perfect chameleon effect, even for bleached teeth
- Perfect gloss due to optimal light reflection
- Permanent colour fidelity
- Highly resistant to discolouration



10112 | OMNICHROMA  
Syringe (à 4 g)

10122 | OMNICHROMA  
20 capsules (à 0.2 g)

10113 | OMNICHROMA BLOCKER  
Syringe (à 4 g)


10123 | OMNICHROMA BLOCKER  
20 capsules (à 0.2 g)

1 Composite  
2 Viscosities  
All Benefits



10232 | OMNICHROMA FLOW  
Syringe (à 3 g)

10234 | OMNICHROMA BLOCKER FLOW  
Syringe (à 3 g)



TOKUYAMA Dental Deutschland GmbH  
Am Landwehrbach 5 • 48341 Altenberge  
Germany

Tel: +49 2505 938513 • Fax: +49 2505 938515  
[info@tokuyama-dental.de](mailto:info@tokuyama-dental.de)  
[www.omnichroma.eu](http://www.omnichroma.eu)

From the Technology Pioneer.