





These are our standard textiles and colors. For more information, please visit our website at [vela-medical.com](http://vela-medical.com).


## Illusion




BLUE 68746 (VELA NO. 063)



GREY 57841 (VELA NO. 059)



BLACK 00110 (VELA NO. 058)




LIGHT GREY 57069 (VELA NO. 067)

### Illusion


**Cleaning:** Wipe regularly with warm water and neutral detergent. Use a clean white cloth or a sponge.

**Disinfection:** Use detergents containing up to 75% ethanol or chlorine disinfection containing active chlorines up to 5000 ppm. After cleaning, rinse with cold, clean water for longer durability.




MINK GREY 46508 (VELA NO. 062)

## Dolce



BLACK - 6728 (VELA NO. 771)




METEOR GREY - 6727 (VELA NO. 777)

### Dolce

**Cleaning:** Wipe with a solution of 10% mild soap and warm water using a soft cloth. Rinse with a clean, damp cloth.

**Disinfection:** Use any of the following solutions: ethanol, water and bleach (5:5:2) or bleach and water (1:5) or ethanol and water (1:1).



PEACOCK BLUE - 6723 (VELA NO. 770)