











Overview: ICX-Prosthetics Box Instruments








Article-Nr.		Article	Material	Price €
960001		ICX-Ratchet, adjustable torque of 10-45Ncm	Stainless steel	159*
950098		ICX-Hex tool 1.4 long, Length 28mm for all abutment screws, healing and cover caps	Stainless steel	29*
C-015-100025		ICX-Hex tool 1.4 medium, for all abutment screws	Stainless steel	49*
960002		Hex tool 0.9 for t-ecco female inner part and ICX-mini cover screw	Stainless steel	29*
C-015-100030		ICX-Hexagon instrument (one piece) SW 1.2mm , medium, L 27mm	Stainless steel	49*
C-015-100007		ICX-t-bona insertion tool (ISO-shaft) Installation of all t-bona screws	Stainless steel	49*
C-020-000001		ICX-multi insertion tool SW 3.2mm for straight multi abutments	Stainless steel	49*
960007		ICX-Universal manual adapter to combine all instruments with elbow receptacle (ISO shaft), for the torque ratchet	Stainless steel	49*
960004		ICX-Universal drill extension for all ICX drills usable up to max. 35Ncm	Stainless steel	49*
C-015-100005		ICX-Extraction tool, L 33mm, for stripping tight ICX abutments	Stainless steel	49*

ICX-Shop: www.medentis.de






ICX
The FAIR Implant-System

Overview: ICX Accessories

ICX and ICX-ACTIVE-MASTER Insertion Tools












Article-Nr.		Article	Material	Price €
C-015-100006		ICX insertion tool with ISO-shaft with 1mm T-Bona screw exclusively (standard). For all ICX- <i>mini</i> implants Ø 2.9mm	Stainless steel	49,-*
C-015-100001		ICX insertion tool, short, length 23mm Insertion tool short with ISO-shaft For all Ø 3.45/3.75/4.1/4.8 implants	Stainless steel	29,-*
C-015-100002		ICX insertion tool, medium, length 29mm Insertion tool medium sized with ISO-shaft For all Ø 3.45/3.75/4.1/4.8 implants	Stainless steel	29,-*
C-015-100004		ICX insertion tool long, length 40mm Insertion tool long with with ISO-shaft For all Ø 3.45/3.75/4.1/4.8 implants	Stainless steel	29,-*
C-015-100021		ICX insertion tool, medium, Insertion tool medium sized with ISO-shaft For all Ø 3.45/3.75/4.1/4.8 implants	Stainless steel	49,-*
C-015-100040		ICX insertion tool with ISO-shaft Lockable, length 30mm The innovative click system allows a simple and safe removal of the ICX implants.	Stainless steel	59,-*
C-015-100035		ICX insertion tool, adjustable The innovative click system enables a simple and safe removal of the ICX-implants.	Stainless steel	79,-*

ICX – Other Accessories

Article-Nr.		Article	Material	Price €
C-015-100009		Countertool for <i>ICX-mini</i> (for ball head removal)	Stainless steel	29*
C-014-230315		ICX parallel pin Ø 2.3/Ø 3.15 for drilling cavity control	Titanium	29*
C-014-350410		ICX parallel pin Ø 3.5/Ø 4.1 for drilling cavity control	Titanium	29*
C-014-006002		ICX bone profiler for the manual removal of crestal bone that may grow over the implant platform during healing	Stainless steel	9.90*
C-015-100020		ICX surgical driver, for instruments with ISO-shaft	Stainless steel	99*

Overview: ICX Accessories

ICX-Hex Tools

Article-Nr.		Article	Material	Price €
950099		ICX-Hex tool 1.4 medium , length 24mm for all abutment screws	Stainless steel	29*
950098		ICX-Hex tool 1.4 long , length 28mm for all abutment screws	Stainless steel	29*
950097		ICX-Hex tool 1.4 extra long , length 45mm, for all abutment screws	Stainless steel	29*
C-015-100023		ICX-Hex tool SW 1.4 short , one piece complete	Stainless steel	49*
C-015-100025		ICX-Hex tool one piece SW 1.4 medium , for all abutment screws	Stainless steel	49*
C-015-100024		ICX-Hex tool SW 1.4 long , piece complete	Stainless steel	49*
C-015-100033		ICX-Hex tool SW 1.2 ISO-shaft	Stainless steel	29*
C-015-100030		ICX-Hex tool SW 1.2 medium piece complete	Stainless steel	49*
C-015-100036		ICX-Hex tool one piece SW 1.2mm , L 32mm	Stainless steel	49*
960002		Hex tool 0.9 for t-ecco female inner part and ICX-mini cover screw	Stainless steel	29*
C-015-100029		ICX-Hex tool one piece SW 0.9mm , length 38.5mm for t-ecco abutments and ICX-mini cover screw	Stainless steel	49*

ICX
The FAIR Implant-System